


GalvoSurge Dental Implant Cleaning System GS 1000 Implants compatibility list

The GalvoSurge Dental Implant Cleaning System is suitable for all electroconductive dental implants (titanium) independent from their geometries, surface or manufacturer if the implant shoulder corresponds to the following dimensions:



Below you can find a list of implants within the above stated dimensions that can be cleaned with the GalvoSurge Dental Implant Cleaning System GS 1000. This list is **non-extensive**, GalvoSurge can be used on many more implant systems if they are within the dimensions indicated above.

Implant Company	Implant Line	Implant Platform	Implant Diameter
Biomet 3i®	Certain Internal	Certain Internal 3.4	3.25
		Certain Internal 4.1	4.0
		Certain Internal 5.0	5.0
	External Hex		3.25
			3.75
			4.0
Camlog®	Camlog K	Camlog K 3.3	3.3
		Camlog K 3.8	3.8
		Camlog K 4.3	4.3
		Camlog K 5.0	5.0
		Camlog K 6.0	6.0
	Camlog J	Camlog J 3.3	3.3
		Camlog J 3.8	3.8
		Camlog J 4.3	4.3
		Camlog J 5.0	5.0
		Camlog J 6.0	6.0
	Conelog	Conelog 3.3	3.3
		Conelog 3.8	3.8
		Conelog 4.3	4.3
		Conelog 5.0	5.0

Implant Company	Implant Line	Implant Platform	Implant Diameter
Camlog®	Isy	Isy	3.8
			4.4
			5.0
Dentsply Astra Tech Implant System™	"Classical" AstraTech TX Line	3.0 (XS)	3.0
		3.5/4.0 (S)	3.5
			4.0
		4.5/5.0 (L)	4.5
		5.0	
	EV Line	3.0	3.0
		3.6	3.6
		4.2	4.2
4.8		4.8	
Dentsply Ankylos®	Older C line without index	C Type	3.5
			4.5
			5.5
			7.0
	Younger X line with index	X Type	3.5
			4.5
			5.5
			7.0
Dentsply Sirona XiVE®	S Line	3.0 (Friadent)	3.0
		3.4 (Friadent)	3.4
		3.8 (Friadent)	3.8
		4.5 (Friadent)	4.5
		5.5 (Friadent)	5.5
Dentsply Sirona Frialit®	Frialit	3.4 (Friadent)	3.4
		3.8 (Friadent)	3.8
		4.5 (Friadent)	4.5
		5.5 (Friadent)	5.5
Nobel Biocare	NobelActive®, NobelParallel™ CC, NobelReplace® CC, NobelReplace® CC PMC	Conical Connection - 3.0	3.0
		Conical Connection - NP	3.5
			3.75
		Conical Connection - RP	4.3
			5.0
	Conical Connection - WP	5.5	
	NobelReplace® Tapered, Replace Select™ Tapered, Replace Select™ Tapered PMC	Tri-Channel Connection - NP	3.5
		Tri-Channel Connection - RP	4.0
			4.3
		Tri-Channel Connection - WP	5.0
6.0			

Implant Company	Implant Line	Implant Platform	Implant Diameter
Nobel Biocare	Branemark System®, NobelSpeedy® Groovy	External Hex Connection - NP	3.3
		External Hex Connection - RP	3.75
			4.0
		External Hex Connection - WP	5.0
			6.0
SIC® invent	Ace Line	3.3	3.4
		3.3	4.0
		4.2	4.5
		4.2	5.0
Straumann®	Tissue Level	NN	3.5
		NNC	3.3
		RN	3.3
			4.1
			4.8
	WN	4.8	
	Bone Level	NC	3.3
		RC	4.1
4.8			

GMT70879 © Nobel Biocare Services AG, 2020. All rights reserved. Distributed by: Nobel Biocare. Legal manufacturer: GalvoSurge Dental AG, Nöllenstrasse 15a, CH-9443 Widnau, Switzerland. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to nobelbiocare.com/trademarks for more information. Product images are not necessarily to scale. All product images are for illustration purposes only and may not be an exact representation of the product. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed clinician, medical professional or physician. See Instructions For Use for full prescribing information, including indications, contraindications, warnings and precautions.