

# Paths to NobelProcera® restorations

## Intra-oral scanners

Vendor	Model	Supported NobelProcera® restorations			Software or program used			Scan locator used		
		Cement retained	Screw-retained	Implant bars	DTX Studio	Scan & Design	Open Access	Elos Medtech	Nobel Biocare	3Shape
Carestream	3500	Single-unit 2-6 units	Single-unit		●	●	●	✓		
	3600	Single-unit 2-6 units	Single-unit 2-3 units		●	●	●	✓		
iTero	I/II Generation	Single-unit 2-6 units	Single-unit		●	●	●	✓		
	Elements	Single-unit 2-6 units	Single-unit 2-3 units		●		●	✓		
Medit	i500 i700	Single-unit 2-14 units	Single-unit 2-3 units		●	●	●	✓		
Planmeca	Emerald S	Single-unit 2-6 units			●	●	●	✓		
Sirona	CEREC Omnicam	Single-unit 2-6 units	Single-unit		●	●	●	✓		
	Primescan	Single-unit 2-14 units	Single-unit 2-3 units		●		●	✓		
	Bluecam Apollo DI	Single-unit 2-6 units	Single-unit		●	●	●	✓		
3M	True Definition True def mobile version	Single-unit 2-6 units	Single-unit		●	●	●	✓		
3Shape	Trios Mono/Color Trios 2/3/4	Single-unit 2-6 units	Single-unit 2-3 units		●	●	●	✓		✓

# Desktop scanners

Vendor	Model	Supported NobelProcera® restorations			Software or program used			Scan locator used		
		Cement retained	Screw-retained	Implant bars	DTX Studio	Scan & Design	Open Access	Elos Medtech	Nobel Biocare	3Shape
Nobel Biocare	LS 3	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	
	2G	Single-unit 2-14 units	Single-unit 2-14 units	✓	●		●		✓	
3Shape	D640	Single-unit	Single-unit		●	●	●	✓	✓	✓
	D500, D700, D710/720/750, D800/810/850, D900/L, D1000/2000, E1, E2, E3, E4	Single-unit 2-14 units	Single-unit 2-14 units	✓	●			✓	✓	
		Single-unit 2-14 units	Single-unit 2-14 units		●	●	●	✓	✓	✓
Amann Girrbach	Map 200	Single-unit 2-14 units	Single-unit 2-14 units		● <sup>1</sup>	●	●	✓	✓	
	Map 400/600	Single-unit 2-14 units	Single-unit 2-14 units	✓	● <sup>1</sup>	●	●	✓	✓	
Dentalwings	7 series	Single-unit 2-14 units	Single-unit 2-14 units		●	●		✓	✓	
DOF	Freedom HD	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	
Medit	Identica T300/T500, t500	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	
Sirona	inEos X5	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	
SmartOptics	Activity 888	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	
Zirkonzahn	S600	Single-unit 2-14 units	Single-unit 2-14 units	✓	●	●	●	✓	✓	

<sup>1</sup> Supported by DTX Studio Lab 1.10 release or newer



[nobelbiocare.com/lab](https://www.nobelbiocare.com/lab)

GMT 78156 GB 2201 © Nobel Biocare Services AG, 2022. All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare. Please refer to [nobelbiocare.com/trademarks](https://www.nobelbiocare.com/trademarks) for more information. Product images are not necessarily to scale. Disclaimer: Some products may not be regulatory cleared/released for sale in all markets. Please contact the local Nobel Biocare sales office for current product assortment and availability. For prescription use only. Caution: Federal (United States) law restricts this device to sale by or on the order of a licensed dentist. See Instructions for Use for full prescribing information, including indications, contraindications, warnings and precautions.

